

AHN myTip RB

Rack Box System without Tips



Effortless and secure handling - easy one hand opening and closing of rack box / secure rack box closing / rack boxes can be stacked for easy storage

Universal fit - racks are compatible with multi channel pipettes (standard 9 mm SBS footprint)

Sterile handling - racks and rack boxes are autoclavable (121 °C / 20 min)

Rack boxes available in small and large size / racks in 8x12 format and 6x10 format



metal free

CE IVD

autoclavable
121°C
20 min

Ordering information

Volume	Format	Colour	Sterilized	Pack Type	Sales Unit	Cat. No.
10/20 µL	8x12	clear		carton	5x10	4-500-C5-0
200 µL	8x12	clear		carton	5x10	4-501-C5-0
200/300 µL	8x12	clear		carton	5x10	4-502-C5-0
1000 µL, long	8x12	clear		carton	5x10	4-503-C5-0
1000 µL, short	6x10	clear		carton	5x10	4-504-C5-0



+49(0)3631/46594-04
+49(0)3631/46594-10
info@ahn-bio.de

Uthleber Weg 14
D-99734 Nordhausen
www.ahn-bio.de

ahn®
biotechnologie GmbH